

# **Polymer-based pre-filled syringes designed to minimize the aggregation risk of sensitive biodrugs**

Nelson Labs Virtual Symposium  
Small Volume Parenteral Packaging  
30-31 March 2022




# OVERVIEW

- ✓ Corporate Introduction
- ✓ Challenges of biopharmaceuticals
- ✓ Polymer-based pre-fillable syringes as a system approach for sensitive biodrugs
- ✓ Silicone oil-free solution
- ✓ Minimizing the risks of protein oxidation
- ✓ Aspects of contact materials: extractables & leachables

# Introduction



# Terumo was founded by medical scientists



© Courtesy of the Kitasato Institute archives

Dr. Shibasaburo Kitasato  
(1853 – 1931)


- Terumo was founded by several medical scientists led by Dr. Shibasaburo Kitasato.
- Discovered the immune antibody against the tetanus toxin and established a serum therapy for tetanus.
- Identified a plague bacillus and paved the way for preventive medicine.



<b>Head Office:</b>	Shibuya-ku, Tokyo, Japan
<b>Foundation:</b>	September, 1921
<b>Stock:</b>	1st Section of Tokyo Stock Exchange
<b>CapEx:</b>	0.72 billion US\$
<b>Net Sales:</b>	5.75 billion US\$
<b>Employees:</b>	26 400




## Terumo's Sales by region



# Terumo's Business

## Sales by Business



\*Approx \$5.75 bn (1 USD =106.754 JPY)

### Cardiac and Vascular Company

- Vascular Intervention
- Intervention oncology
- Cardiovascular surgery



### General Hospital Company

- Hospital Systems Division
- Alliance Business  
(Collaborating with  
pharmaceutical manufacturers)



### Blood Management Company

- Blood centers
- Therapeutic apheresis and cell collection
- Cell processing



# Challenges of biopharmaceuticals




# Challenges of sensitive biodrugs

## Defects in quality can cause fatal incidents

- Peginesatide (Omontys®; Affymax, Inc.,) was voluntarily withdrawn from the market less than a year after launch.
- 49 cases of anaphylaxis, including 7 fatalities, were reported.
- Sub-visible particles linked to the hypersensitivity cases.

**Joseph Kotarek, Christine Stuart, Silvia H. De Paoli, et al.** Subvisible Particle Content, Formulation, and Dose of an Erythropoietin Peptide Mimetic Product Are Associated With Severe Adverse Postmarketing Events. *Journal of Pharmaceutical Science*, Nov 23 2005






## Immunogenicity is a concern

Various factors that affect immunogenicity are indicated in FDA Guidance

- **Protein aggregates** may elicit or enhance immune responses
- **Chemical modifications**, such as oxidation may elicit immune responses
- **Sub-visible particles** in the size range 0.1-10 microns have a strong potential to be immunogenic
- **Leached materials** from the container closure system may be a source of materials that enhance immunogenicity by chemically modifying the therapeutic protein product




Guidance for Industry; Immunogenicity  
Assessment for Therapeutic Protein  
Products (August 2014)

# New chapters and revisions in USP

Revisions were made in response to further availability of biopharmaceuticals

- US pharmacopeia (USP) new chapter
  - **USP <1663>** Extractables good practice
  - **USP <1664>** Assessment of leachables
- USP revisions
  - **USP <661.1> and <661.2>**  
Related to plastic packaging materials
- USP new informational chapter
  - **USP <1787>**  
Measurement of sub-visible particulate matter in therapeutic protein injections, intended to supplement USP <787>



# Factors potentially affecting biopharmaceutical PFS quality

**Various causing factors = need a system approach**

Factor	Issue	Caused by
Physical stress	Aggregation by silicone oil	Silicone oil
	Aggregation by tungsten	Manufacturing process
	Aggregation by glue	Manufacturing process
	Aggregation by excess physical stress	Head space
Chemical stress	Denaturation by leachables from container closure	Leachables
	Oxidation by dissolved oxygen	Dissolved oxygen
	Oxidation by remaining radicals	Irradiation
Others	Breakage	Component
	Delamination	Component
	Particles	Manufacturing process

# Polymer PFS as a system approach



# A system approach with polymer-based pre-fillable syringes

## Designed for sensitive biodrugs

- Polymer barrel (COP: Cyclo Olefin Polymer)
  - Insert molding → **Tungsten-free, Glue-free**
  - High break resistance
- Prevent protein aggregation
  - i-coating™ technology → **Silicone oil-free**
  - Low sub-visible particles
- Prevent protein oxidation
  - Secondary packaging with oxygen absorber
  - Steam sterilization → **No radical generation**
- Low extractable syringe system
  - COP, Chlorinated butyl rubber (CIIR)
  - Steam sterilization → Low extractable






# Silicone oil-free system: Stopper coating technology

i-coating™ on the plunger stopper achieves the silicone oil-free system





i-coating™ provides a smooth surface on PLA-JEX™ stopper that works in combination with the dimensional tolerance of the polymer barrel to provide an enhanced container closure integrity.

Uncoated stopper




i-coating™ stopper




# i-coating™ stopper

Silicone oil-free system shows low sub-visible particle count



# Oxygen permeability

Higher permeability of polymer can be an advantage




The value of polymer permeation and combined use of oxygen absorber with the packaging system  
(a proven concept with >15 years of experience)

# Packaging structure for elimination of dissolved oxygen


Prevent protein oxidation by dissolved oxygen

The de-oxygenated package system



# Effect of de-oxygenated packaging system


## Enables elimination of the dissolved oxygen





# Advantage of steam sterilization

Prevent protein oxidation by radicals




# Low EXTRACTABLE syringe system

## Achieved by combining material selection and applying steam sterilization

### STUDY 1

Extraction condition: PLAJEX filled with WFI at 121°C for 60 min  
Analysis: Organic extractable compounds by LC/MS

#### Material comparison




# Low LEACHABLE syringe system

## STUDY 2

Leachable study condition: 5ml PLA JEX filled with WFI at 25°C / 60% RH for 36 months

PP label with acrylic-based glue


Analysis: Organic leachable compounds & elements



- Detected levels of organic leachables were low and those components were not considered as a concern of toxicity.
- Label-related leachables could not be detected.

# Low LEACHABLE syringe system

## STUDY 2



- Detected levels of elements leachables were low and those components were not considered as a concern of toxicity.
- ✓ No differential compounds were found compared with the control
- ✓ There were no leachables detected related to the secondary packaging and label.
- ✓ Leachable levels from the primary packaging were much lower than toxicity threshold.

# Low LEACHABLE syringe system

## STUDY 3

Leachable study condition: 1ml PLA JEX filled with Humira simulant at 2-8°C / darkness for 36 months

PP label with acrylic-based glue

Analysis: Organic leachable compounds, elements & acrylic acid

Analytical type	Compound	Toxicity Class	(µg/day)	Origin
VOC	Ethyl acetate	50 mg/day <sup>1)</sup>	2.24	Unconfirmed
SVOC	Cyclohexanone	50 µg/day <sup>2)</sup>	0.19	Unconfirmed
	2-Ethyl-1-hexanol	50 µg/day	0.29	Unconfirmed
	Pentadecanoic acid, dimethyl ester	50 µg/day	0.09	Unconfirmed
	2-(2-Butoxyethoxy)ethanol	50 µg/day	0.08	Unconfirmed
	Dimethyl adipate	50 µg/day	0.07	Unconfirmed
	C <sub>21</sub> H <sub>40</sub> Rubber oligomer	50 µg/day	0.08	Rubber (Tip cap or Plunger stopper)
	4-Undecylbenzene sulfonic acid	50 µg/day	2.08	Rubber (Tip cap or Plunger stopper)
NVOC	4-Dodecylbenzene sulfonic acid	50 µg/day	4.16	
	4-Tridecylbenzene sulfonic acid	50 µg/day	7.52	
	BHT (3,5-ditert-butyl-4-hydroxytoluene)	50 µg/day	0.05	
Elements	Fe	-	0.16	All
	Mg	-	0.62	All
	Mn	-	0.02	Rubber (Tip cap or Plunger stopper)
	Si	-	1.12	Rubber (Tip cap or Plunger stopper)

\*1) ICH Q3C guideline, \*2) PQRI guideline

- Leachable levels from the PLA JEX container closure system with Humira simulant were much lower than the toxicity threshold.
- Label-related leachables could not be detected.
- ✓ PLA JEX offers a low leachable PFS system in applying low extractable / leachable COP and ClIR materials and well-selected secondary packaging components.



# Conclusion

## PLAJEX™ with i-coating™ stopper

- Risks for sensitive biodrugs such as aggregation and oxidation can be minimized by a system approach with polymer-based pre-fillable syringes.

- Tungsten-free, Glue-free
- High break resistance
- Silicone oil-free
- Low sub-visible particles
- Secondary packaging with oxygen absorber
- Steam sterilization = Radical-free
- Low extractable materials and sterilization method



# References

1. **Functional Evaluation and Characterization of a Newly Developed Silicone Oil-Free Prefillable Syringe System**, KEISUKE YOSHINO *et al*, Wiley Online Library online, 18 March 2014
  2. **Assessment of the effects of sterilization methods on protein drug stability by elucidating decomposition mechanism and material analysis**, Koji Nakamura *et al*, International Journal of Pharmaceutics, 2014
- Impact of Sterilization Method on Protein Aggregation and Particle Formation in Polymer-Based Syringes**, Hideaki Kiminami *et al*, Journal of Pharmaceutical Sciences, 2016

