

The Necessity of Extractable and Leachable Qualifications for Lyophilized Drug **Products: Some Fallacies Addressed**

Sona Kovackova, PhD **E&L Expert SVP, Nelson Labs Europe**



Table of Content

- 1. Lyo Containers and Administration Devices: what are we talking about?
- 2. Lyo Primary Containers: What do the US & EU Guidelines tell us?
- 3. The Interaction Mechanism between the Lyo Drug Product and its Primary Packaging: **OUTGASSING**
- 4. The Challenges and Consequences of Leachables Testing after Reconstitution
- 5. Reactive Leachables a Concern for LYO-DP Applications
- 6. Key Learnings



1. Lyo Containers and Administration Devices: what are we taking about?





1. Lyo Containers and Administration Devices: what are we taking about?



Drug Administration Sets

(Not discussed in this presentation)







2. Lyo Primary Containers: What do the US & EU Guidelines tell us?





Lyo Primary Containers: What do the US & EU Guidelines tell us?



Revision of "Table 1" in USP <1664>, Originally Included into the FDA Guidance for Industry (1999): "Container/Closure systems for Packaging Human Drugs and Biologics"









The "Low" Likelihood of Packaging - DP Interaction for LYO SVP's: based upon the observation that:

- 1. the *interactio*n between a solid (Lyo cake) a material (eg rubber) *is limited*
- 2. AND, there is *limited direct contact* between Lyo cake and Rubber closure
- However the Interaction Mechanism between LYO Cake and Components may not need direct contact (see later)
- BE CAREFUL when "rationalizing" a LYO application as being Non Critical!!!





Lyo Primary Containers: What do the US & EU Guidelines tell us? The EM(E)A Guideline on "Plastic Immediate Packaging



For solid active substances and solid dosage forms: the risk of interaction is low and generally does not require a content/container interaction study. Solid dosage forms intended for inhalation or parenteral use, e.g. lyophilised products, may need interaction studies between the packaging material and the components of the formulation.





Lyo Primary Containers: What do the US & EU Guidelines tell us?

Conclusion wrt Regulatory Guidance (US, EU)

POTENTIAL AMBIGUITY in the **NEED** and/or subsequent **DEPTH** of E/L evaluations for Lyo-Containers







3. The Interaction Mechanism between the Lyo Drug Product and its Primary Packaging: OUTGASSING





The Interaction Mechanism: During Long Term STORAGE of the LYO-Cake



© 2022 Nelson Labs NV | All Rights Reserved.



Outgassing of rubbers is mainly an issue for:



Rubber Stopper



Outgassing of rubbers is mainly an issue for: Volatile Organic Compounds

LYO CAKE: EXTREMELY GOOD ADSORBENT

ACCUMULATION OF VOC/SVOC LEACHABLES **Potential Interactions** (Reactive Leachables)

ADSORPTION ON LYO-CAKE SURFACE



The Interaction Mechanism: During RECONSTITUTION of the Lyophilized DP

RECONSTITUTION





of adsorbed leachables in reconstitution solution

Mainly Volatile & Semi-Volatile Organic compounds (OUTGASSING)



Reconstituted (Liquid) DP Liquid Interaction with **Primary Packaging Components**

Although SHORT TERM Contact, also Non-Volatile Compounds, Metals and lons may be released





Basic Analytical Techniques



The Extractable Study

Good Identification Practices for **Organic Extractables and Leachables Via Mass Spectrometry**

All Rights Reserved. © 2022 Nelson Labs NV

The Leachable Study

© 2022 Nelson Labs NV | All Rights Reserved.

SCREENING FOR UNEXPECTED LEACHABLES

EXAMPLE OF FULL LEACHABLE STUDY

Type of Solution	Storage Time (Months)				
	0	3	6	12	24
Drug Product in Rubber Sealed Vials (Test Item) at 5 \pm 3 $^\circ$ C	×	×	×	×	×
Drug Product in Inert Containers (Blank) at 5 \pm 3 $^{\circ}$ C	×	×	×	×	×
Drug Product in Rubber Sealed Vials (Test Item) at 25 \pm 3 $^\circ$ C	_	×	×	-	-
Drug Product in Inert Containers (Blank) at 25 \pm 3 $^{\circ}$ C	-	×	×	-	-
× = sampling time point		5			

The Interaction Mechanism: During ADMINISTRATION of the RECONSTITUTED DP

TRANSIENT TO SHORT TERM CONTACT WITH DRUG ADMINISTRATION DEVICES

Transient or Short Term: No Long Term Accumulation However: Potentially Increased Levels of VOC, SVOC, NVOC, Metals and Ions Potentially: Include the Administration Step in the Leachable Procedures

IS IT A MEDICAL DEVICE, OR IS IT A COMBINATION PRODUCT HOWEVER, THESE CONSIDERATIONS WILL NOT BE ADDRESSED IN THIS PRESENTATION

4. The Challenges and Consequences of Leachables Profiling after Reconstitution

The Challenges and Consequences of Leachables Profiling after Reconstitution

Blank solution: What defines a good blank and what is the issue with LYO?

A good blank solution is a **leachables free drug product**!!

A good blank solution is best from the **same drug product batch** as the contact samples

A good blank solution is **put on controlled storage together with the contact samples**

20

EXAMPLE: SCREENING LEACHABLE WITHOUT BLANK REFERENCE = IMPURITIES PROFILING

No Blank for LYO-product is available: "you know what? I will screen without the Blank Baseline"

75 Organic Impurities = More than only Leachables from the Primary Packaging!

© 2022 Nelson Labs NV All Rights Reserved.

21

The Challenges and Consequences of Leachables Profiling after Reconstitution

THE CHALLENGE:

IN LYO-APPLICATIONS: **NO TRUE AGED BLANK DRUG PRODUCT IS AVAILABLE**

POTENTIAL SOLUTIONS:

THERE IS NO "PERFECT" SOLUTION Consider t = 0 as blank baseline throughout the complete study Disadvantage:

- Some Compounds may already be present after a short term contact
- Stability / Degradation of DP is not accounted for

Use upright position as a "reference" for inverted position Disadvantage:

• *Migration Mechanism* is not based upon material contact, rather base release of VOC/S-VOC into headspace of vial with subsequent adsorption

Use "Lab Prepared & Freshly Prepd DP", containing all the ingredient Disadvantage:

- Some volatile impurities may be present in "Prepared DP" but not in material (lyophilization also removes a lot of VOC/S-VOC impurities)
- Degradation of DP may/will be different in solution compared to Lyoph

ANOPHIN Emptying Lyo-vial right after manufacturing and store the dry LYO-pro-Disadvantage:

- Find the right container for long term storage?
- How to keep the LYO-DP dry over the shelf life?

PREFERRE

ATAL

NOOPTON

ANDPHON.

The Challenges and Consequences of Leachables Profiling after Reconstitution

LEACHABLE STUDY DESIGN

RELY MORE ON TARGET ANALYSIS OF LEACHABLES IN THE DRUG PRODUCT WITH QUANTITATIVE METHODS

© 2022 Nelson Labs NV All Rights Reserved.

EVALUATE REACTIVITY

WITH SMALL MOLECULE API'S

Adduct Formation of an API with the C₁₃H₂₃Br and C₂₁H₃₉Br oligomers

© 2022 Nelson Labs NV | All Rights Reserved.

26

WITH EXCIPIENTS AN OTHER INGREDIENTS

Lyophilized Drug Product in a glass vial with a rubber stopper (without coating):

Excipients a.o.: Glycine

WITH EXCIPIENTS AN OTHER INGREDIENTS

Lyophilized Drug Product in a glass vial with a rubber stopper (without coating):

Excipients a.o.: Polysorbate 20

C21H39Br rubber oligomer source leachable from a rubber stopper

© 2022 Nelson Labs NV | All Rights Reserved.

Fatty acids source –

Schiff base reaction 29

WITH THERAPEUTIC PEPTIDES AND PROTEINS

See presentation of Paulo Forte on Reactive Leachables and the R&D work we did and why it is relevant for Imunogenicity assessment for Therapeutic Protein Nelson Labs.

A Sotera Health company

© 2022 Nelson Labs NV | All Rights Reserved.

Experimental R&D work, performed by Nelson Labs on Lyophilized Glargine

6. Key Learnings

Key Learnings

- The (regulatory) need or depth of performing an E/L-qualification is ambiguous
- The interaction mechanism between the LYO DP and the rubber stopper is **NOT** based upon a direct contact mechanism but goes through the mechanism of OUTGASSING
- This may lead to accumulation of leachables, adsorbed onto the Lyo cake
- In addition adsorbed Leachables may be reactive and may lead to DP degradation
- Adding a screening step in your formal leachable studies can evaluate the occurrence of leachable reactivity
- This can assist in mitigating the risk of the occurrence of immunogenicity through chemical modification of the **Therapeutic Protein**
- Consider the right rubber quality for your Lyophilized C/C-system consult supplier

Thank you

Questions? InfoEurope@nelsonlabs.com +32 16 40 04 84

Register for FREE access

for this presentation and much more expert content on

Soterahealth.com/academy

Safeguarding Global Health®

Expert Lab Testing & Advisory Services

nelsonlabs.com sales@nelsonlabs.com +1 801-232-6293

Comprehensive Sterilization Solutions & Expert Advisory Services

sterigenics.com +1 800-472-4508

Cobalt-60

nordion.com service@nordion.com +1 800-465-3666

© 2022 Nelson Labs NV | All Rights Reserved.

Nordion.

A Sotera Health company

Reliable Global Supply of