

Chemical Interactions between Leachables and Biopharmaceuticals

COMBINED IN-SILICO AND EXPERIMENTAL MODEL (WITH INSULIN AS A MARKER COMPOUND)
TO MONITOR THE POTENTIAL IMPACT ON QUALITY AND SAFETY OF THERAPEUTIC PROTEINS

Paulo Forte
Study Director
Extractables & Leachables



Eprex case




- Eprex = Human Recombinant EPO
- Introduced in late '80 early '90 – Janssen Cilag
- Increase Red Blood Cell-count in Chronic Kidney Disease Patients
- Until '98: no side effects
- From '98 onwards: increased incidence of PRCA
 - *Caused a drop in Hematocrit (instead of an increase)*
 - *Immune response*

Eprex case: Incidents of AB-mediated Pure Red Cell Aplasia

CONTAINER
CLOSURE
SYSTEM

Antibody-mediated
PRCA



Red Blood Cell levels substantially reduced because
of an Anti-body mediated immune response


Peak ^a	Compound	Average concentration ^b
1	Unknown	Unknown
2	Bisphenol A	0.070
3	4- <i>tert</i> -amylphenol	0.046
4	2-chloro-4- <i>tert</i> -amylphenol	0.037
5	Vultac [®] 2 disulfide	0.778
6	2,2'-methylene-bis-4- <i>tert</i> -amylphenol	0.243
7	Vultac [®] 2 trisulfide	0.235
8	Vultac [®] 2 tetrasulfide	0.142
9	Vultac [®] 2 pentasulfide	0.063
10	Vultac [®] 2 hexasulfide	0.024

QUESTION: Who could have predicted an anti-body mediated immune response, based upon those analytical data?

Mode of action - Hypothesis in the early work:


Leachables (one or more) could have **adjuvant-like properties**, “boosting” an immune response, which is **causing ADA’s (Anti-Drug-Antibodies) to be formed**

ADA’s attacked both endogenous & exogenous EPO
ultimately resulting in **a substantial decrease of Red Blood Cells (PRCA/Anemia)**



However, the “**adjuvant like properties**” of the detected compounds were studied in **animal models**, but **no ADA’s were observed**.

Mode of action – New Line of Thinking:



**Reactive Leachables may form covalent bonds with
Biologics and this may lead to Immuno Responses**

Relevance of reactive leachables

8. *Container Closure Considerations*

Leached materials from the container closure system may be a source of materials that enhance immunogenicity, either by chemically modifying the therapeutic protein product or by having direct immune adjuvant activity, including the following:

- Organic compounds with immunomodulatory activity may be eluted from container closure materials by polysorbate-containing formulations: a leachable organic compound involved in vulcanization was found in a polysorbate formulated product when the stopper surfaces were not Teflon coated (Boven et al. 2005).

—————> Refers to Eprex case

Guidance for Industry

Immunogenicity Assessment for
Therapeutic Protein Products

U.S. Department of Health and Human Services
Food and Drug Administration
Center for Drug Evaluation and Research (CDER)
Center for Biologics Evaluation and Research (CBER)

August 2014
Clinical/Medical

Risk evaluation

Recommendations

Whenever possible, sponsors should obtain detailed information regarding a description of all raw materials used in the manufacture of the container closure systems for their products. Sponsors should conduct a comprehensive extractables and leachables laboratory assessment using multiple analytical techniques to assess the attributes of the container-closure system that could interact with and degrade protein therapeutic products.

Testing for leachables should be performed on the product under stress conditions,⁹ as well as under real-time storage conditions, because in some cases the amount of leachables increases dramatically over time and at elevated temperatures. Product compatibility testing should be performed to assess the effects of container closure system materials and all leachables on product quality.

Guidance for Industry

Immunogenicity Assessment for Therapeutic Protein Products

U.S. Department of Health and Human Services
Food and Drug Administration
Center for Drug Evaluation and Research (CDER)
Center for Biologics Evaluation and Research (CBER)

August 2014
Clinical Medical




... Conduct E/L assessment with multiple techniques to assess the attributes of the C/C system that **could interact with/degrade the protein therapeutic product...**



Immunogenicity risk evaluation in E&L – how?


Immunogenicity risk

Time consuming





Protein (family)
dependent

Immunogenicity risk evaluation in E&L – how?




Initial Toxicological
evaluation

Immunogenicity risk evaluation in E&L – how?




QUESTION 1:

CAN THE REACTIVITY OF LEACHABLES
WITH PROTEINS AND PEPTIDES BE PREDICTED BY AN
IN SILICO MODEL?




Can we predict reactive leachables?




Leached materials from the container closure system may be a source of materials that enhance immunogenicity, either by chemically modifying the therapeutic protein product or by having direct immune adjuvant activity, including the following:

(FDA guidance 2014)


in silico – reactivity (models based on)



$C_{13}H_{23}Br$ – rubber oligomer



acrylic acid



palmitic acid

glutathione

potency by GSH

extremely reactive

highly reactive

moderately reactive

slightly reactive

suspect

lysine

DPRA (13%) Lys

reactive

grey zone

unreactive

cysteine

DPRA (13%) Cys

reactive

grey zone

unreactive

NELSON LABS

high reactivity


moderate reactivity

low reactivity


in silico prediction

QSAR TOOLBOX


in silico – reactivity (models based on)



$C_{13}H_{23}Br$ – rubber oligomer



acrylic acid



palmitic acid

glutathione

potency by GSH

extremely reactive

highly reactive

moderately reactive

slightly reactive

suspect

lysine

DPRA (13%) Lys

reactive

grey zone

unreactive

cysteine

DPRA (13%) Cys

reactive

grey zone

unreactive

NELSON LABS

high reactivity

moderate reactivity

low reactivity

in silico prediction


QSAR TOOLBOX

CAS	Name	Nelson reactivity
57-10-5	palmitic acid	
84-66-2	diethyl phthalate	
5323-65-9	2-chloro-4-tert-amylphenol	
	vultac 3	
96-76-4	2-4-di-tert-butylphenol	
78-42-2	tris(2-ethylhexyl) phosphat	
93-89-0	ethyl benzoate	low
95-47-6	o-xylene	low
96-37-7	methylcyc	low
78-59-1	alpl	moderate
80-56-8		moderate
544-10-5		moderate
1120-72		moderate
7		moderate
	acid	high
	er oligomer	high
	etophenone	high
	horone diisocyanate	high
9	styrene oxide	high
85	phthalic anhydride	high
80	methyl methacrylate	high

**AT NELSON LABS:
REACTIVITY ASSESSMENT OF OVER 3500
COMPOUNDS**

QUESTION 2:


CAN THE PREDICTED REACTIVITY OF
LEACHABLE COMPOUNDS BE OBSERVED IN
BIOLOGIC DRUG PRODUCTS?



Why Lyophilized Insulin (Glargine)?

- **High relevancy** because of widespread use and increasing delivery device market
- **Lyophilized** Insulin is **readily available** as reference material (EDQM)
- **Acceptable analytical complexity**


Wet lab experiments: leachables



Leachable	In silico prediction?	Reacted?	Leachable	In silico prediction?	Reacted?
Palmitic acid	Low reactivity	?	<i>cis</i> -9,10-epoxystearic acid	High reactivity	?
Diethyl phthalate	Low reactivity		C ₁₃ H ₂₃ Br rubber oligomer	High reactivity	
2-chloro-4-tert-amylphenol	Low reactivity		2-Chloroacetophenone	High reactivity	
Vultac 3	Low reactivity		Isophorone diisocyanate	High reactivity	
Caprolactam	Low reactivity		2,6-di-tert-butyl-p-benzoquinone	High reactivity	
7,9-di-tert-butyl-1-oxaspiro [4.5]deca-6,9-diene-2,8-dione	Low reactivity		Diocetyl disulfide	High reactivity	


Wet lab experiments: overview

Insulin Glargine + leachables



Which parts react?


Structure Insulin Glargine



How to analyze?

Instrumentation

Q-Exactive Orbitrap Liquid Chromatography Mass Spectrometry




How to prepare?

Sample preparation

Lyophilized Insulin is spiked with leachables
Solvent dried off, incubation 4 days @ 40°C

Insulin Glargine Blank reference solution




11jan2019_001 #2002-2016 RT: 17.99-18.11 AV: 15 SB: 2 17.66, 18.56 NL: 6.92E8
T: FTMS + p ESI Full ms [150.0000-2000.0000]

m/z 1010.812, 6+
Molecular mass 6063.863 Da

6+


High resolution
Accurate mass






22



Insulin Glargine Adduct



Insulin Glargine + leachables (2)



Wet lab experiments: result overview



Leachable	In silico prediction?	Reacted?	Leachable	In silico prediction?	Reacted?
Palmitic acid	Low reactivity	✗	<i>cis</i> -9,10-epoxystearic acid	High reactivity	✓
Diethyl phthalate	Low reactivity	✗	C ₁₃ H ₂₃ Br rubber oligomer	High reactivity	✓
2-chloro-4-tert-amylphenol	Low reactivity	✗	2-Chloroacetophenone	High reactivity	✓
Vultac 3	Low reactivity	✗	Isophorone diisocyanate	High reactivity	✓
Caprolactam	Low reactivity	✗	2,6-di-tert-butyl-p-benzoquinone	High reactivity	✓
7,9-di-tert-butyl-1-oxaspiro [4.5]deca-6,9-diene-2,8-dione	Low reactivity	✓	Diocetyl disulfide	High reactivity	✗

*Oxaspiro: low reactivity predicted in OECD toolbox
However: OECD toolbox also indicated a Michael Addition as a potential reaction mechanism, which was observed through a low reactivity*

*Diocetyl Disulfide: HIGH reactivity predicted in OECD toolbox
However: OECD tests were performed in a reducing environment
Creating sulfide radicals that are very reactive.
Our testing was not performed under such conditions...*

How to include in a leachable study?



extractable study
find potential leachables



initial risk
assessment



leachables study



PREDICT TOXICITY

NELSON's FIT SCREEN
(Fast Initial Toxicological QSAR Screening)

supports material choice

indication of potential
toxicological risk to patient – **direct risk**



PREDICT REACTIVITY

Additional

supports material choice



get compounds for stability study – **direct risk**

indication of potential
immunogenic risk to patient – **indirect risk**



For Medical Device E/L: this tool could be
used as a predictor for sensitizing properties

***In silico* prediction of reactivity:**

- cost effective with over 3500 chemical structures assessed on their reactivity
- Understanding the risk of reactive leachables can assist in
 - Optimizing Stability Studies
 - Risk mitigation for Immunogenicity
 - Material Selection

CONTAINER CLOSURE SYSTEM
COMPATIBLE?



THANK YOU

I would like to thank the following people who all were involved in the scientific work & interpretation:

- Dr. Piet Christiaens
- Dr. Ank Reumer
- Ruud Cuyvers
- Dr. Philippe Verlinde
- Dr. Lise Vanderkelen

Thank you

Questions?

InfoEurope@nelsonlabs.com

+32 16 40 04 84



SoteraHealth



ACADEMY

Register for **FREE** access

for this presentation
and much more expert content on

Soterahealth.com/academy



Sotera[®]Health

Safeguarding Global Health[®]



Nelson Labs[®]

A Sotera Health company

**Expert Lab Testing &
Advisory Services**

nelsonlabs.com
sales@nelsonlabs.com
+1 801-232-6293



Sterigenics[®]

A Sotera Health company

**Comprehensive
Sterilization Solutions &
Expert Advisory Services**

sterigenics.com
+1 800-472-4508



Nordion[®]

A Sotera Health company

**Reliable Global Supply of
Cobalt-60**

nordion.com
service@nordion.com
+1 800-465-3666