



USP Water Test Matrix

	Particulate Matter	Bacterial Endotoxins	Sterility	Microbial Limits	Antimicrobial Agents	Conductivity	Total Organic Carbon	Oxidizable Substances	pH, Sulfate, Calcium, Carbon Dioxide	Ammonia, Chloride
Purified Water						X	X			
Sterile Purified Water			X			X		X		
Pure Steam		X				X	X			
Water for Injection		X				X	X			
Bacteriostatic Water for Injection	X	X	X		X				X	
Sterile Water for Injection	X	X	X					X	X	X
Sterile Water for Inhalation		X	X			X		X		
Sterile Water for Irrigation		X	X			X		X		
Water for Hemodialysis		X		X		X		X		