

Nelson Laboratories, LLC
Process Challenge Device Specification Sheet

This information is specific o PCD orders related to lab services.

Please complete a PCD Order Form and Submit to GP-PCDOrders@nelsonlabs.com

PCD Type	PCD Description	BI Type(s)	BI manufacturer(s)	Estimated D-value
NLI-PCD-001	Consists of a Tygon tubing insert inside a syringe with a vented cap, packaged in a single Tyvek [®] pouch.	BI Strip	Steris NAMSA Mesa Labs SPS Medical	18.8
NLI-PCD-002	Consists of a syringe with a non-vented cap packaged in a single Tyvek [®] pouch	BI Strip Self-Contained BI	Steris NAMSA Mesa Labs SPS Medical 3M	36.8
NLI-PCD-003	Consists of a syringe with vented cap packaged in a single Tyvek [®] pouch	BI Strip Self-Contained BI	Steris NAMSA Mesa Labs SPS Medical 3M	3.3
NLI-PCD-004	Consists of a syringe with a non-vented cap packaged in a single Tyvek [®] pouch	BI Strip Self-Contained BI	Steris NAMSA Mesa Labs SPS Medical 3M	62.4
NLI-PCD-005	Consists of a BI packaged in a single Tyvek [®] pouch	BI Strip Self-Contained BI	Steris NAMSA Mesa Labs SPS Medical 3M	2.8
NLI-PCD-006	Consist of a BI inside a double Tyvek [®] pouch	BI Strip Self-Contained BI	Steris NAMSA Mesa Labs SPS Medical 3M	3.1
NLI-PCD-007	Consists of Tygon tubing with vinyl caps packaged in a single Tyvek [®] pouch	BI Strip	Steris NAMSA Mesa Labs SPS Medical	20.2
NLI-PCD-008	Consists of Tygon tubing with vinyl caps packaged in a single Tyvek [®] pouch	BI Strip	Steris NAMSA Mesa Labs SPS Medical	14.9

NLI-PCD-009	Consists of a syringe with no cap packaged in a single Tyvek [®] pouch	BI Strip Self-Contained BI	Steris NAMSA Mesa Labs SPS Medical 3M	2.5
NLI-PCD-010	Consists of a syringe with a Tygon tube and vinyl cap attached with a luer lock packaged in a single Tyvek [®] pouch	BI Strip Self-Contained BI	Steris NAMSA Mesa Labs SPS Medical 3M	27.2
NLI-PCD-011	Consists of a syringe with a Tygon tube and vinyl cap attached with a luer lock packaged in a single Tyvek [®] pouch	BI Strip Self-Contained BI	Steris NAMSA Mesa Labs SPS Medical 3M	14.9
NLI-PCD-012	Consists of a syringe with a non-vented cap packaged in a single Tyvek [®] pouch	BI Strip	Steris NAMSA Mesa Labs SPS Medical	35.6
NLI-PCD-013	Consists of a syringe with a non-vented cap packaged in a single Tyvek [®] pouch	BI Strip Self-Contained BI	Steris NAMSA Mesa Labs SPS Medical 3M	42.3