



FaceMask Testing Parameters Request:

1. Testing Parameters

Which side to test (if sending in flat sheets of material):

Inside Outside

2. Synthetic Blood Penetration

Please let us know what test pressure you would like to use:

N/A 80 mm HG 120 mm HG 160 mm HG

3. Particle Filtration Efficiency

Which particle size are we to use:

N/A 0.1 microns 0.3 microns 0.5 microns 1.0 microns

Note: The particle size required by ASTM F2100, Standard Specification for Performance of Materials Used in Medical Face Masks is 0.1 microns.

*Include this form with shipment **along with** the sample submission forms.